



DiaSorin Italia S.p.A.
Via Crescentino snc - 13040 Saluggia (VC) - Italy
www.diasorin.com
Tel. +39.0161.4871



Changes: §4, §6;
Deletions: -

LIAISON® Treponema Screen (REF 310840)

1. INTENDED USE

The LIAISON® Treponema Screen assay uses chemiluminescent immunoassay (CLIA) technology for the qualitative determination of specific total antibodies to *Treponema pallidum* in human serum or plasma samples. The test has to be performed on the LIAISON® Analyzer family.

2. SUMMARY AND EXPLANATION OF THE TEST

Syphilis is a disease, usually sexually transmitted, caused by infection with the spirochaetae bacterium *Treponema pallidum*. Congenital transmission of *Treponema pallidum* by transplacental passage from infected mothers and contagion through blood transfusion have been also reported in the literature. Infection is systemic from the outset and the disease is characterized by periods of latency, often in excess of twenty years. Syphilis natural course is divided into three phases. After an incubation period lasting about three weeks, a non-painful skin lesion (chancre) appears, often associated with regional lymphadenopathy (primary phase). The disease progresses into a secondary, disseminated, phase, accompanied by general mucocutaneous lesions and lymphadenopathy. If a *Treponema pallidum* infection is allowed to progress into the late phases of the disease, the secondary phase is followed by a period of subclinical infection (latent syphilis) that can be detected only by serological tests, and by a late or tertiary phase, observed only in a small number of patients, characterized by progressive disease.

Syphilis may be diagnosed by different serological laboratory tests, which are also useful in establishing the stage of the disease when adopted together with other clinical tests. Serological diagnosis is generally established by using two standard tests combined, non-treponemal antibody tests for screening purposes and specific treponemal antibody tests for confirmation purposes. The two most commonly used non-treponemal antibody tests are the Venereal Disease Research Laboratory test (VDRL) and the Rapid Plasma Reagin test (RPR), the latter a simplified variation of VDRL.

In many laboratories, the screening technique of choice is the *Treponema pallidum* haemagglutination test (TPHA) that consists of agglutination of erythrocytes coated with specific *Treponema pallidum* antigens and detects specific IgG and IgM antibodies. However, interpretation of TPHA results is subjective and cannot be completely automated. Different immunoenzymatic kits employing bacterial lysate and, later on, purified recombinant proteins have been recently placed on the market.

Positive results can be confirmed by a reference test, the fluorescence treponemal antibody absorption test (FTA-Abs), that allows detection of both IgG and IgM. Its performance, however, is laborious and time-consuming, unsuitable for screening large numbers of samples, and interpretation of results is subjective.

In contrast to non-treponemal tests, antigen-specific treponemal test reactivity persists after treatment, often for life. Blood and plasma donors are screened for *Treponema pallidum* antibodies using both non-treponemal and treponemal tests.

3. PRINCIPLE OF THE PROCEDURE

The method for determination of specific total antibodies to *Treponema pallidum* is a one-step sandwich chemiluminescence immunoassay (CLIA). Recombinant antigens specific for *Treponema pallidum* are used for coating magnetic particles (solid phase) and are linked to an isoluminol derivative (isoluminol-antigen conjugate). During the incubation, *Treponema pallidum* antibodies present in calibrators, samples or controls bind to the solid phase as well as to the antigen conjugate. After the incubation, the unbound material is removed with a wash cycle.

Subsequently, the starter reagents are added and a flash chemiluminescence reaction is thus induced. The light signal, and hence the amount of isoluminol-antigen conjugate, is measured by a photomultiplier as relative light units (RLU) and is indicative of *Treponema pallidum* total antibody concentration present in calibrators, samples or controls.

4. MATERIALS PROVIDED

Reagent integral

Magnetic particles (2.3 mL)	[SORB]	Magnetic particles coated with <i>Treponema pallidum</i> recombinant antigens (obtained in <i>E. coli</i>), BSA, PBS buffer, < 0.1% sodium azide.
Calibrator 1 (1.4 mL)	[CAL1]	Human serum/plasma containing low <i>Treponema pallidum</i> antibody levels, BSA, phosphate buffer, 0.2% ProClin™ 300, an inert yellow dye. The calibrator concentrations are referenced to an in-house antibody preparation.
Calibrator 2 (1.4 mL)	[CAL2]	Human serum/plasma containing high <i>Treponema pallidum</i> antibody levels, BSA, phosphate buffer, 0.2% ProClin™ 300, an inert blue dye. The calibrator concentrations are referenced to an in-house antibody preparation.
Specimen diluent (13 mL)	[DILSPE]	Proteins, EDTA, phosphate buffer, 0.2% ProClin™ 300, an inert blue dye.
Conjugate (9 mL)	[CONJ]	<i>Treponema pallidum</i> recombinant antigens (obtained in <i>E. coli</i>), conjugated to an isoluminol derivative, BSA, PBS buffer, 0.2% ProClin™ 300, preservatives.
Number of tests		200

All reagents are supplied ready to use. The order of reagents reflects the layout of containers in the reagent integral.

Materials required but not provided (system related)

LIAISON® XL Analyzer	LIAISON® Analyzer
LIAISON® XL Cuvettes ([REF] X0016). LIAISON® XL Disposable Tips ([REF] X0015). LIAISON® XL Starter Kit ([REF] 319200). LIAISON® Wash/System Liquid ([REF] 319100). LIAISON® XL Waste Bags ([REF] X0025). –	LIAISON® Module ([REF] 319130). – LIAISON® Starter Kit ([REF] 319102) or LIAISON® XL Starter Kit ([REF] 319200). LIAISON® Light Check 12 ([REF] 319150). LIAISON® Wash/System Liquid ([REF] 319100). LIAISON® Waste Bags ([REF] 450003). LIAISON® Cleaning Kit ([REF] 310990).

Additionally required materials

LIAISON® Treponema Screen controls (negative and positive) ([REF] 310841).

5. WARNINGS AND PRECAUTIONS

For *in vitro* diagnostic use.

All serum and plasma units used to produce the components provided in this kit have been tested for the presence of HBsAg, anti-HCV, anti-HIV-1, anti-HIV-2 and found to be non-reactive. As, however, no test method can offer absolute assurance that pathogens are absent, all specimens of human origin should be considered potentially infectious and handled with care.

6. SAFETY PRECAUTIONS

Do not eat, drink, smoke or apply cosmetics in the assay laboratory.

Do not pipette by mouth.

Avoid direct contact with potentially infected material by wearing laboratory clothing, protective goggles, and disposable gloves. Wash hands thoroughly at the end of each assay.

Avoid splashing or forming an aerosol. All drops of biological reagent must be removed with a sodium hypochlorite solution with 0.5% active chlorine, and the means used must be treated as infected waste.

All samples and reagents containing biological materials used for the assay must be considered as potentially able to transmit infectious agents. The waste must be handled with care and disposed of in compliance with the laboratory guidelines and the statutory provisions in force in each Country. Any materials for reuse must be appropriately sterilized in compliance with the local laws and guidelines. Check the effectiveness of the sterilization/decontamination cycle.

Do not use kits or components beyond the expiration date given on the label.

Pursuant to EC Regulation 1272/2008 (CLP) hazardous reagents are classified and labelled as follow:

REAGENTS:	[CAL]1, [CAL]2, [CONJ], [DIL]SPE
CLASSIFICATION	Skin sens. 1A H317 Aquatic chronic 3 H412
SIGNAL WORD:	Warning
SYMBOLS / PICTOGRAMS:	 GHS07 Exclamation mark
HAZARD STATEMENTS:	H317 May cause an allergic skin reaction. H412 Harmful to aquatic life with long lasting effects.
PRECAUTIONARY STATEMENTS:	P261 Avoid breathing dust/fume/gas/mist vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P273 Avoid release to the environment. P362 Take off contaminated clothing and wash before reuse.
CONTAINS: (only substances prescribed pursuant to Article 18 of EC Regulation 1272/2008).	reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1) (ProClin™ 300).

Pursuant to EC Regulation 1272/2008 (CLP), [SORB] is labeled as EUH210 safety data sheets available on request.

For additional information see Safety Data Sheets available on www.diasorin.com.

7. PREPARATION OF REAGENT INTEGRAL

Please note the following important reagent handling precautions:

Resuspension of magnetic particles

Magnetic particles must be completely resuspended before the integral is placed on the instrument. Follow the steps below to ensure complete suspension:

Before the seal is removed, rotate the small wheel at the magnetic particle compartment until the colour of the suspension has changed to brown. Gentle and careful side-to-side mixing may assist in the suspension of the magnetic particles (avoid foam formation). Visually check the bottom of the magnetic particle vial to confirm that all settled magnetic particles have resuspended. Carefully wipe the surface of each septum to remove residual liquid. Repeat as necessary until the magnetic particles are completely resuspended.

An incomplete magnetic particles resuspension may cause variable and inaccurate analytical results.

Foaming of reagents

In order to ensure optimal performance of the integral, foaming of reagents should be avoided. Adhere to the recommendation below to prevent this occurrence:

Visually inspect the reagents, calibrators in particular (position two and three following the magnetic particle vial), to ensure there is no foaming present before using the integral. If foam is present after resuspension of the magnetic particles, place the integral on the instrument and allow the foam to dissipate. The integral is ready to use once the foam has dissipated and the integral has remained onboard and mixing.

Loading of integral into the reagent area

LIAISON® Analyzer

- Place the integral into the reagent area of the analyzer with the bar code label facing left and let it stand for 30 minutes before using. The analyzer automatically stirs and completely resuspends the magnetic particles.
- Follow the analyzer operator's manual to load the specimens and start the run.

LIAISON® XL Analyzer

- LIAISON® XL Analyzer is equipped with a built-in magnetic device which aids in the dispersal of microparticles prior to placement of a reagent integral into the reagent area of the analyzer. Refer to the analyzer operator's manual for details.
 - a. Insert the reagent integral into the dedicated slot.
 - b. Allow the reagent integral to remain in the magnetic device for at least 30 seconds (up to several minutes). Repeat as necessary.
- Place the integral into the reagent area of the analyzer with the label facing left and let it stand for 15 minutes before using. The analyzer automatically stirs and completely resuspends the magnetic particles.
- Follow the analyzer operator's manual to load the specimens and start the run.

8. STORAGE AND STABILITY OF REAGENT INTEGRAL

- **Sealed:** Stable at 2-8°C until the expiry date.
- **Opened on board or at 2-8°C:** Minimum stability four weeks.

After this period, it is still possible to keep on using the reagent integral provided that the controls are found within the expected ranges.
- Use always the same analyzer for a reagent integral already opened.
- Use storage rack provided with the analyzer for upright storage of reagent integral.
- Do not freeze.
- Keep upright for storage to facilitate later proper resuspension of magnetic particles.
- Keep away from direct light.

9. SPECIMEN COLLECTION AND PREPARATION

Either human serum or plasma may be used. The anticoagulants citrate, EDTA, heparin, have been tested and may be used with this assay. Post-mortem specimens, collected up to 24 hours after death, have been tested and may be also used in the assay. The correct specimen type must be used in the assay.

Follow tube manufacturers' instructions carefully when using collection containers. Blood should be collected aseptically by venipuncture and the serum or plasma separated from clot, red cells or gel separator, after centrifugation.

Centrifugation conditions range from 1,000 to 3,000 g for 10 minutes. Conditions may vary depending on tube manufacturers recommendations. Use of alternate centrifugation conditions should be evaluated and validated by the laboratory.

Before shipping specimens, serum or plasma specimens should be removed from clot, red cells or gel separator. Specimens may be shipped in dry ice (frozen), in wet ice (for 2°-8°C) or at room temperature (20°-25°C), by following sample storage limitations described below.

Uncontrolled transport conditions (in terms of temperature and time) can cause inaccurate analytical results.

During validation studies, specimen collection tubes commercially available at the time of testing were used. Therefore not all collection tubes from all manufacturers have been evaluated. Blood collection devices from various manufacturers may contain substances which could affect the test results in some cases (Bowen et al., Clinical Biochemistry, 43, 4-25, 2010).

Concerning storage limitations, if the assay is performed within seven days of sample collection, the samples removed from red cells, clot or gel separator may be kept at 2°-8°C; otherwise they should be aliquoted and stored deep-frozen (-20°C or below). Eight samples with different reactivity were stored for seven days at 2°-8°C and underwent four freeze-thaw cycles. The results showed no significant differences; however multiple freeze-thaw cycles should be avoided. If samples are stored frozen, mix thawed samples well before testing.

Samples removed from red cells, clot or gel separator having particulate matter, fibrin, turbidity, lipaemia, or erythrocyte debris, specimens that have been stored at room temperature (20°-25°C), or frozen and thawed, or samples requiring repeat testing, require clarification by further centrifugation (it's recommended 10,000 g for 10') before testing, to improve consistency of results. Specimens with a lipid layer on the top should be transferred in a secondary tube, taking care to transfer only the clarified material. Grossly haemolyzed or lipaemic samples as well as samples containing particulate matter or exhibiting obvious microbial contamination should not be tested. Check for and remove air bubbles before assaying.

The minimum volume required is 230 µL specimen (80 µL specimen + 150 µL dead volume).

10. CALIBRATION

Test of assay specific calibrators allows the detected relative light unit (RLU) values to adjust the assigned master curve. Each calibration solution allows four calibrations to be performed.

Recalibration in triplicate is mandatory whenever at least one of the following conditions occurs:

The analyzer should be calibrated in triplicate whenever one of the following conditions occurs:

- A new lot of reagent integral or of Starter Kit is used.
- The previous calibration was performed more than two weeks before.
- The analyzer has been serviced.
- Control values lie outside the expected ranges.

LIAISON® Analyzer: Calibrator values are stored in the bar codes on the integral label.

LIAISON® XL Analyzer: Calibrator values are stored in the Radio Frequency IDentification transponder (RFID Tag).

11. ASSAY PROCEDURE

Strict adherence to the analyzer operator's manual ensures proper assay performance.

LIAISON® Analyzer. Each test parameter is identified via the bar codes on the reagent integral label. In the event that the barcode label cannot be read by the analyzer, the integral cannot be used. Do not discard the reagent integral; contact your local DiaSorin technical support for instruction.

LIAISON® XL Analyzer. Each test parameter is identified via information encoded in the reagent integral Radio Frequency IDentification transponder (RFID Tag). In the event that the RFID Tag cannot be read by the analyzer, the integral cannot be used. Do not discard the reagent integral; contact your local DiaSorin technical support for instruction.

The analyzer operations are as follows:

1. Dispense specimen diluent, coated magnetic particles and conjugate into the reaction module.
2. Dispense calibrators, controls or specimens.
3. Incubate.
4. Wash with Wash/System liquid.
5. Add the Starter reagents 1 and 2 and measure the light emitted.

12. QUALITY CONTROL

LIAISON® controls should be run in singlicate to monitor the assay performance. Quality control must be performed by running LIAISON® Treponema Screen controls

- (a) at least once per day of use,
- (b) whenever a new reagent integral is used,
- (c) whenever the kit is calibrated,
- (d) whenever a new lot of Starter Reagents is used,
- (e) to assess adequacy of performance of the open integral beyond four weeks, or in agreement with guidelines or requirements of local regulations or accredited organizations.

Control values must lie within the expected ranges: whenever one or both controls lie outside the expected ranges, calibration should be repeated and controls retested. If control values obtained after successful calibration lie repeatedly outside the predefined ranges, the test should be repeated using an unopened control vial. If control values lie outside the expected ranges, patient results must not be reported.

The performance of other controls should be evaluated for compatibility with this assay before they are used. Appropriate value ranges should then be established for quality control materials used.

13. INTERPRETATION OF RESULTS

The analyzer automatically calculates *Treponema pallidum* total antibody levels expressed as index value and grades the results. For details, refer to the analyzer operator's manual.

Calibrators and controls may give different RLU or dose results on LIAISON® and LIAISON® XL, but patient results are equivalent.

The cut-off value discriminating between the presence and the absence of *Treponema pallidum* total antibodies has an index value of 1. Sample results should be interpreted as follows:

Samples with *Treponema pallidum* antibody levels below an index value of 0.9 should be graded *negative*.

Samples with *Treponema pallidum* antibody levels ranging between an index value of 0.9 and 1.1 should be graded *equivocal*.

Samples with *Treponema pallidum* antibody levels equal to or above an index value of 1.1 should be graded *initially reactive*.

A sample resulting *initially reactive* or *equivocal* at the first assay should be assayed in duplicate in order to confirm the initial result.

If a sample results repeatedly reactive in at least one replicate it should be considered positive.

Samples which are non-reactive at the second test should be considered negative.

A second sample should be collected and tested no less than one week later when the result is repeatedly equivocal.

A negative result for total antibodies to *Treponema pallidum* generally indicates that the patient has not been infected, but does not always rule out acute syphilis, because the infection may be in its very early stage and the patient may be still unable to synthesize *Treponema pallidum* specific antibodies, or the antibodies may be present in undetectable levels. If clinical exposure to *Treponema pallidum* is suspected despite a negative or equivocal finding, a second sample should be collected and tested later during the course of infection.

A positive result for total antibodies to *Treponema pallidum* generally indicates exposure to the pathogen (acute or past infection). A single specimen, however, can only help estimate the serological status of the individual.

14. LIMITATIONS OF THE PROCEDURE

This test allows screening for the presence of *Treponema pallidum* total antibodies. It detects both recent and past infections, but it cannot distinguish between different antibody classes.

Detection of *Treponema pallidum* total antibodies may indicate recent, past or successfully treated syphilis: this test, therefore, cannot discriminate between active and treated disease and, as a consequence, may not be used for determining response to therapy.

In order to monitor the effectiveness of therapy clinical and serologic evaluation (non-treponemal tests) should be performed.

LIAISON® Treponema Screen test may produce positive results that score negative with non-treponemal tests (VDRL, RPR), because it detects *Treponema pallidum* antibodies that persist for life. RPR tests usually produce negative results in past infections, because they detect heterophilic antibodies that are present only in the early phase of infection.

A non reactive test result does not exclude the possibility of exposure to or infection with syphilis. *T. pallidum* antibodies may be undetectable in some stages of the infection and in some clinical conditions.

A skillful technique and strict adherence to the instructions are necessary to obtain reliable results.

Bacterial contamination or heat inactivation of the specimens may affect the test results.

Test results are reported qualitatively as positive or negative for the presence of *Treponema pallidum* total antibodies. However, diagnosis of infectious diseases should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures as well as in association with medical judgement. Therefore a precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. Test results, however, must be interpreted with caution in immunocompromized individuals since their antibody levels may be affected by this condition.

Integrals may not be exchanged between analyzer types (LIAISON® and LIAISON® XL). Once an integral has been introduced to a particular analyzer type, it must always be used on that analyzer until it has been exhausted. Due to traceability issues resulting from the above statement, patient follow-ups may not be concluded between analyzer types. These must be accomplished on one particular analyzer type (either LIAISON® or LIAISON® XL).

Before testing cadaveric specimens, collection and centrifugation procedures should be carefully applied. After death, haemolysis and other changes (including proteolysis and dilution) occur in blood, which may lead to False Negative and False Positive in testing. In subjects transfused immediately prior to death high percentage of haemodilution can affect the performance of the test due to analyte dilution.

15. SPECIFIC PERFORMANCE CHARACTERISTICS

15.1. Analytical specificity

Analytical specificity may be defined as the ability of the assay to accurately detect specific analyte in the presence of potentially interfering factors in the sample matrix (e.g., anticoagulants, haemolysis, effects of sample treatment), or cross-reactive antibodies.

Interference. Controlled studies of potentially interfering substances or conditions showed that the assay performance was not affected by anticoagulants (sodium citrate, EDTA, heparin), haemolysis (up to 1000 mg/dL haemoglobin), lipaemia (up to 3000 mg/dL triglycerides), bilirubinaemia (up to 20 mg/dL bilirubin), or by freeze-thaw cycles of samples.

Cross-reactions. As a rule, the presence of potentially cross-reactive antibodies does not interfere in the assay. The antibodies investigated were: (a) immunoglobulins to various infectious agents – such as hCMV, EBV, VZV, rubella virus, *Borrelia burgdorferi*, *Treponema denticola* or *Toxoplasma gondii* – (b) anti-nuclear (ANA) antibodies and rheumatoid factor (anti-Fc immunoglobulin) antibodies.

15.2. Precision with LIAISON® Analyzer

Different samples, containing different concentrations of specific analyte, were assayed to estimate repeatability and reproducibility of the assay (i.e., within- and between-assay variability). The variability shown in the tables below did not result in sample misclassification. **The results refer to the groups of samples investigated and are not guaranteed specification as differences may exist between laboratories and locations.**

Repeatability. Twenty replicates were performed in the same run to evaluate repeatability.

Repeatability	A	B	C	D
Number of determinations	20	20	20	20
Mean (index value)	0.13	8.17	51.48	58.53
Standard deviation	0.02	0.22	1.39	1.63
Coefficient of variation (%)	13.1	2.6	2.7	2.8

Reproducibility. Several replicates were performed in different days (one or two runs per day) with three different lots of integral to evaluate reproducibility. The tests were performed in two sites, in house and in an independent laboratory.

Reproducibility	SITE 1				SITE 2			
	E	F	G	H	E	F	G	H
LOT No. 01								
Number of determinations	20	20	20	20	11	11	11	11
Mean (index value)	0.18	11.02	40.20	56.45	0.15	12.50	41.68	57.47
Min value	0.14	9.40	26.90	47.00	0.10	11.30	38.10	42.40
Max value	0.20	12.80	46.70	61.00	0.19	14.00	47.00	66.70
Coefficient of variation (%)	10.8	7.3	12.1	5.9	15.4	7.0	7.1	11.2
LOT No. 02								
Number of determinations	20	20	20	20	11	11	11	11
Mean (index value)	0.13	12.24	38.08	56.20	0.07	13.15	40.00	57.36
Min value	0.09	11.00	33.20	49.00	0.05	11.90	36.70	50.70
Max value	0.26	13.90	44.30	62.30	0.09	14.00	43.00	64.00
Coefficient of variation (%)	39.2	6.3	8.5	5.9	21.9	4.6	7.9	6.3
LOT No. 03								
Number of determinations	20	20	20	20	11	11	11	11
Mean (index value)	0.06	12.60	40.83	54.45	0.03	12.57	38.85	55.80
Min value	0.04	11.40	34.70	47.30	0.03	11.10	33.70	49.70
Max value	0.10	13.50	49.30	61.30	0.03	13.60	42.30	62.00
Coefficient of variation (%)	24.7	3.7	12.1	8.4	0.0	5.7	9.1	6.6

15.3. Precision with LIAISON® XL Analyzer

Different samples, containing different concentrations of specific analyte, were assayed to determine repeatability and reproducibility of the assay (i.e., within- and between-assay variability). The variability shown in the tables below did not result in sample misclassification. **The results refer to the groups of samples investigated and are not guaranteed specification as differences may exist between laboratories and locations.**

Repeatability. Twenty replicates were performed in the same run to evaluate repeatability.

Repeatability	1	2	3	Positive control
Number of determinations	20	20	20	20
Mean (index value)	0.565	3.07	3.99	11.1
Standard deviation	0.051	0.056	0.082	0.34
Coefficient of variation (%)	9.1	1.8	2.0	3.0
Min. value	0.492	2.99	3.80	10.5
Max. value	0.656	3.19	4.12	11.6

Reproducibility. Twenty replicates were performed in different days (one or two runs per day) to evaluate reproducibility.

Reproducibility	1	2	3	Positive control
Number of determinations	20	20	20	20
Mean (index value)	0.328	2.91	3.78	10.5
Standard deviation	0.053	0.10	0.12	0.47
Coefficient of variation (%)	16.0	3.6	3.2	4.5
Min. value	0.261	2.77	3.56	9.89
Max. value	0.422	3.17	4.06	11.2

15.4. High-dose hook effect

Whenever samples containing extremely high antibody concentrations are tested in a one-step sandwich method, the hook effect can mimic concentrations lower than real.

Analysis of hook effect was evaluated by testing four high-titred samples positive for total antibodies to *Treponema pallidum*. All samples resulted in a concentration value above the assay range that would be expected with high-titred sera, indicating no sample misclassification.

15.5. Diagnostic specificity and sensitivity

Diagnostic specificity and sensitivity were estimated by testing 3690 specimens from different populations. The specimens were tested by several comparison methods and consensus between them as well as the available clinical and serological data were applied to define the expected results. Global results are summarized herebelow; detailed results follow.

Diagnostic specificity: 99.91% (3506/3509) - 95% confidence interval: 99.75-99.98%.

Diagnostic sensitivity: 99.40% (167/168) - 95% confidence interval: 96.73-99.98%.

BLOOD DONORS. Out of 2494 samples from a non-selected blood donor population, 2491 results were negative, one result was equivocal and two results were positive. Diagnostic specificity was 99.92% (2491/2493), calculated upon exclusion of the equivocal result (95% confidence interval: 99.71-99.99%).

CLINICAL SAMPLES. A panel of 131 samples collected from patients in different syphilis phases (7 primary phase, 31 secondary phase, 77 latent syphilis, 5 cardiovascular syphilis, 6 neurosyphilis, 5 congenital syphilis) were characterized by commercially available tests, such as RPR, TPHA, enzyme immunoassay and in-house Western blot, and were evaluated by LIAISON® Treponema Screen test. Six results were equivocal and then excluded, because the comparison tests were not in agreement. Out of 125 classified positive samples, 125 results were positive by LIAISON® Treponema Screen test. Diagnostic sensitivity was 100% (125/125), calculated upon exclusion of equivocal results (95% confidence interval: 97.09-100%).

CROSS-REACTIONS. 65 potentially cross-reactive samples were also tested: out of them, 30 results were positive for *Borrelia burgdorferi* antibodies, 10 results were positive for EBV antibodies, 10 results were positive for *Streptococcus β-haemolyticus* antibodies, 15 results were positive for *Treponema denticola* antibodies. All those potentially cross-reactive samples were negative by LIAISON® Treponema Screen test.

PROSPECTIVE SAMPLES. 1000 samples were collected from non-selected laboratory routine specimens with relevant serological results (commercial enzyme immunoassay and in-house Western blot). As the comparison tests were not in agreement on five results, they were classified equivocal. Considering the results given by enzyme immunoassay and Western blot combined as consensus, diagnostic specificity was 99.89% (950/951 - 95% confidence interval: 99.42-100%) and diagnostic sensitivity was 97.67% (42/43 - 95% confidence interval: 87.71-99.94%), calculated upon exclusion of equivocal results.

15.6. Performance characteristics of cadaveric specimen testing

Performance characteristics of cadaveric specimens testing was determined by testing, according PEI validation protocol*, post-mortem specimens collected up to 24 hours after death in comparison to living donor specimens. 20 post-mortem samples were tested as unspiked and spiked at 2 levels: low positive and medium/high positive. The same procedure was performed with the same number of normal human serum from living donors, tested in parallel as reference to compare with post-mortem sample results. The results obtained were analyzed through calculation of percentage difference between mean of living donors results and mean of post-mortem results, at each reactivity level. In this study, the obtained percentage difference was equal or below 6,9% for each of the tested reactivity levels (see table below). Paired t-test analysis were performed between post-mortem and living donors specimens, spiked at low and medium/high positive levels, demonstrating not significantly difference on two groups (p value <0.05).

Repeatability was assessed using one post-mortem and one living donor specimens, spiked up to a low-level of reactivity with a human serum reactive for antibodies to *Treponema pallidum*. Each specimen was assessed in six replicates in the same run. The obtained percent coefficient of variation (CV%) did not exceed 15%. As reported in the table below 1.5% for the cadaveric specimen and 1.3% for the living donor were found in the study. The results refer to the group of investigated samples and are not guaranteed specifications, as differences may exist between laboratories and locations.

	Sample	Test results	Recovery (%)	t-test	CV%
		Means (index value)	Post-mortem/Living donors	p value	6 replicates
Neat	Post-Mortem unspiked	0.182	n.a.	n.a	n.a
	Living donors unspiked	<0.100			
Low Positive	Post-Mortem spiked	2.41	1.3	0.616	1.5
	Living donors spiked	2.38			1.3
Medium/high Positive	Post-Mortem spiked	5.87	6.9	0.104	n.a
	Living donors spiked	5.49			

* Paul Ehrlich Institute - Proposal for the Validation of Anti-HIV-1/2 or HIV Ag/Ab Combination Assays, Anti-HCV-Assays, HBsAg and Anti-HBc Assays for Use with Cadaveric Samples - 08/05/2014

REFERENCES

S.J. NORRIS
Polypeptides of *Treponema pallidum*: progress toward understanding their structural, functional and immunologic roles.
Microb. Reviews, **57** (3) : 750-779 (1993).

V. SAMBRI et al.
Western immunoblotting with five *Treponema pallidum* recombinant antigens for serologic diagnosis of syphilis.
Clin. Diagn. Lab. Immunol., **8** : 534-539 (2001).

A.E. SINGH, B. ROMANOWSKY
Syphilis: review with emphasis on clinical, epidemiologic and some biologic features.
Clin. Microb. Reviews, **12** (2) : 187-209 (1999).

K. A. WORKOWSKI, G. A. BOLAN
Sexually Transmitted Diseases Treatment Guidelines, 2015.
MMWR, **64** (3) : (2015).

200/007-893, 11 - 2024-03